

Push® MetaGrip®

Provides CMC joint stability for pain-free thumb mobility


The Sophisticated Solution to Thumb CMC Osteoarthritis

- Can be worn during work and household activities as well as golf, tennis, knitting, swimming, etc.
- Imbedded, contoured metal insert offers custom fit
- Provides CMC stability during pinch & grip without impeding other joint movement
- Provides long term durability, resisting abrasion
- Can be worn under a glove
- Can be machine washed on gentle cycle at low temperature
- Won't deform if left in hot environment

Measuring / Sizing / Fitting

The Push MetaGrip must be correctly SIZED and SQUEEZED to work!

- To determine the appropriate size, measure circumference of the hand just below the knuckles, excluding the thumb (see below). This indicates hand size, but not necessarily thumb size. Please take this into consideration when choosing your MetaGrip size. If measurement is near the size break point, try on both sizes to determine best fit.
- Ideal fitting around the thumb is accomplished by trying on the Push MetaGrip. Professional measurement and fitting by a therapist / physician is suggested for optimal results.
- Be sure to follow instructions for squeezing metal insert to fit snugly around your thumb.



Push® MetaGrip® Measurement Chart

Select: RIGHT or LEFT

SIZE	INCHES	CENTIMETERS
0 (X-SM)	5 ⁷ / ₈ - 6 ⁷ / ₈	15 - 17.5
1 (SM)	6 ⁷ / ₈ - 7 ³ / ₄	17.5 - 19.5
2 (MED)	7 ³ / ₄ - 8 ⁷ / ₈	19.5 - 22.5
3 (LG)	8 ⁷ / ₈ - 10 ⁵ / ₈	22.5 - 27



The Push MetaGrip must be correctly SIZED and SQUEEZED to work!