


hot therapy products

TheraTemp® moist heat therapy



- place packs in microwave to heat; no water needed!
- produces moist heat by drawing moisture from the air; retains heat in the therapeutic beads
- elastic belts with hook/loop closure hold wrap in place

universal wrap packs with elastic 3 x 27" belt

| | | | |
|---|---------|-----------|-------|
| A | 11-1282 | 6" x 12" | 27.50 |
| B | 11-1281 | 9" x 12" | 35.00 |
| C | 11-1283 | 12" x 16" | 40.00 |
| D | 11-1288 | 14" x 18" | 50.00 |

specialized contour wrap packs with elastic belt

| | | | |
|---|---------|---|-------|
| E | 11-1280 | cervical, 6" x 24" w/ 3 x 27" belt, 2" x 8" strip | 40.00 |
| F | 11-1285 | lumbar, 9" x 18" w/ 3" x 27" belt | 40.00 |
| G | 11-1286 | eye/sinus, 10" x 3.5" w/ 3" x 14" strap | 27.50 |
| H | 11-1287 | wrist, 5.5" x 12" | 30.00 |

economy heat pads



- 3 heat settings
- washable cover
- add water to the sponge for moist heat; remove sponge for dry heat


dry heat pad

| | | | |
|---------|-------|-----------|-------|
| 11-1130 | small | 12" x 15" | 27.50 |
| 11-1131 | large | 12" x 24" | 37.50 |

moist or dry heat pad

| | | | |
|---------|-------|-----------|-------|
| 11-1132 | small | 12" x 15" | 30.00 |
| 11-1133 | large | 12" x 24" | 37.50 |

moist heat pads - analog and digital



- electric pad provides moist heat without adding water (fomentation)
- analog hand controller adjusts temperature
- digital hand controller adjusts temperature and treatment time
- includes flannel cover

analog pads

| | | | |
|---------|--------|----------------------|-------|
| 11-1120 | king | 13" x 27" | 60.00 |
| 11-1121 | medium | 13" x 13" (pictured) | 50.00 |
| 11-1123 | small | 4" x 14" | 45.00 |

digital pads

| | | | |
|---------|----------|----------------------|-------|
| 11-1135 | king | 26" x 14" | 72.50 |
| 11-1136 | medium | 18" x 14" | 60.00 |
| 11-1137 | mini | 19" x 7" | 50.00 |
| 11-1138 | shoulder | 16" x 12" (pictured) | 72.50 |

Relief Pak® Hot Button® compress





instant treatment



- instant heat therapy without hot water or electricity
- great to use on-the-go to warm up muscles and relieve pain



- activate compress by gently flexing disk within the compress
- instantly heats to approximately 130°F (54°C)
- generates heat for at least 30 minutes
- recharge compress by boiling pad
- pad can be recharged at least 100 times



each

bulk (12 each)

| | | | | | | |
|---|----------|-----------|---------|-------|------------|--------|
| I | small | 3½" x 5½" | 11-1025 | 7.50 | 11-1025-12 | 85.00 |
| J | standard | 5" x 15" | 11-1026 | 12.50 | 11-1026-12 | 142.50 |
| K | neck | 11" x 17" | 11-1027 | 17.50 | 11-1027-12 | 200.00 |
| L | oversize | 10" x 11" | 11-1028 | 15.00 | 11-1028-12 | 170.00 |

hot water bottle



- hospital-grade rubber hot water bottle with screw stopper
- 2-quart capacity

| | | |
|------------|-------|-------|
| 11-1140 | each | 7.50 |
| 11-1140-12 | dozen | 82.50 |