# patient safety sensors and alarms

### seat belt sensors





- buckle release has press to release metal buckle
- E-Z release seat belt is ideal for less active patients
- Early Warning belt has extra 12" of Velcro® to remove before removing belt
- alarm (sold separately) sounds when patient opens buckle
- **A** 59-0100 buckle release seat belt 82.50
- **B** 59-0101 EZ release seat belt 70.00
- C 59-0102 Early Warning belt 80.00 59-0109 12" belt extender 15.00
  - **TR2** alarm
- tamper-resistant: no on/off switch

or nurse call

- connect to sensor pads select 85 or 100 db
- volume or silent alarm tone
- ABS mounting holster with mounting clip
- battery check system
- 4 AA batteries with low battery indicator

59-0200 TR2 alarm 59-0200-10 10 each

#### voice alarm

- plays pre-recorded message (of up to 20 seconds) when alarm is triggered

87.50

59

59

765.00

- use magnetic pull cord or connect
- mounting holster
- 4 AA batteries with low battery indicator

I	-	50 (
	a	
n /	2	

- select 85 or 100 db volume or silent alarn tone
- to sensor pads, seatbelts or nurse call system

muicau	Л
59-0250	voice alarm
59-0250-10	10 each



- thin and sensitive sensor; use with lighter patients
- less compression foam for faster recovery when pressure is released
- attaches with sticky strips on pad
- 6-month warranty
- alarms sold separately

59-0110	chair pad	10 x 15"	45.00
59-0110-10	10 each		405.00
59-0111	bed pad	11 x 30"	85.00
59-0111-10	10 each		765.00

## IQ Duo plus alarm

functions as pull cord and sensor alarm



- silent setting can be used in conjunction with
- select 85 or 100 db volume/alarm tone
- tamper-resistant: no on/off switch; hidden reset button
- mounting holster with strap
- 3 AAA batteries with low battery indicator

-0210	IQ Duo Plus alarm	70.
-0210-10	10 each	630.

## IQ mini pull cord alarm

- magnetic pull cord releases in any direction
- easily attaches to wheelchair or bedside
- non-replaceable battery included

# proximity pad sensors



proximity sensor

- works with patient on chair or bed
- provides earlier fall warning
- works on all air support surfaces
- no weight requirements; one size fits all, cut sensor to desired length
- thin sensor is undetectable to patients
- alarms sold separately

59-0120	Mylar proximity sensor	50.00
59-0120-10	proximity sensor, 10 each	450.00
59-0121	universal connector	67.50
59-0122*	cable (specify suffix)*	35.00

\* specify cable when ordering (030, 030R, 040, 050R, 060R, 070, 110R, 120)

# magnetic pull cord alarm

- magnetic pull cord; pin or connect to nurse call
- select 85 or 100 db volume/alarm tone
- comes with mounting clip
- 2 AAA batteries with >> low battery indicator

59-0230 basic pull cord alarm 37.50 59-0230-10 10 each 360.00

# basic pull pin alarm

- use magnetic cord; pin or connect to nurse call system
- alarm sounds locally and at station when using nurse call feature
- select 85 or 100 db volume/alarm tone
- comes with mounting clip
- 9V battery (not included) with low battery indicator

	voice alarm	105.00	59-0220	mini pull cord alarm	20.00	59-0240	basic pull pin alarm	22.50
)	10 each	900.00	59-0220-10	10 each	180.00	59-0240-10	10 each	202.50

FEI Retail, PO Box 1076, White Plains, New York 10602 USA tel: 888-211-5242 + 914-407-0150 orders@FEIretail.com fax: 888-211-5765 + 914-407-0151



- the nurse call connector







00 .00

