




## gait belts

made of heavy duty 2" wide cotton webbing to provide secure hold for therapists and caregivers while working with patients; belts are adjustable in length

## metal buckle



- nickel chrome plated steel buckle

## quick release plastic buckle



- press top and bottom buttons at the same time to release

## safety quick release buckle



- cover plate prevents accidental opening

50-5130-32	32" belt	9.00	50-5131-32	32" belt	8.00	50-5132-36	36" belt	8.50
50-5130-36	36" belt	9.00	50-5131-36	36" belt	8.00	50-5132-40	40" belt	8.50
50-5130-40	40" belt	10.00	50-5131-40	40" belt	8.50			
50-5130-44	44" belt	10.00	50-5131-44	44" belt	8.50			
50-5130-48	48" belt	10.50	50-5131-48	48" belt	9.00	50-5132-48	48" belt	9.50
50-5130-54	54" belt	10.50	50-5131-54	54" belt	10.00	50-5132-54	54" belt	10.00
50-5130-60	60" belt	11.00	50-5131-60	60" belt	10.50	50-5132-60	60" belt	11.00
50-5130-72	72" belt	11.00	50-5131-72	72" belt	10.50	50-5132-72	72" belt	11.00

## padded transfer belts



- padded transfer belts assist patient transfers
- allows for the comfort and dignity of patient
- pads are 6" wide with four handle loops spaced for safe grip points
- color-coded sizes for quick identification
- available with side release buckle or extra safe auto buckle

side release buckle



auto buckle



## side release buckle


50-5120S	small, blue	24"-30"	40.00
50-5120M	medium, green	28"-48"	42.50
50-5120L	large, black	40"-64"	47.50

## auto buckle

50-5121S	small, blue	24"-30"	50.00
50-5121M	medium, green	28"-48"	52.50
50-5121L	large, black	40"-64"	55.00

## standing / sitting turntable

- rotating pivot disk makes standing or seated transfers easy and safe
- user can sit on turntable and turn on bed or chair
- useful for getting in and out of car seats
- non-skid bottom provides a secure controlled motion
- can be used in conjunction with a gait belt
- 15" diameter, 280 lb capacity



use sitting

use standing



50-5150 turntable 75.00